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## Cross Reactivity

1. The cross reactivity of the 2019-nCoV Antigen Rapid Test Kit was evaluated with a total of 29 microorganisms. None of the microorganisms tested in the following table gave a positive result.

Microorganisms	Concentrations	Microorganisms	Concentrations
Human coronavirus 229E	2.0x10 <sup>5</sup> CID <sub>50</sub> /mL	MERS-coronavirus	2.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL
Human coronavirus OC43	2.0x10 <sup>5</sup> TCID <sub>50</sub> /mL	Chlamydia pneumoniae	2.0 x 10 <sup>6</sup> IFU/mL
Human coronavirus NL63	2.0x10 <sup>5</sup> TCID <sub>50</sub> /mL	Streptococcus pneumoniae	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 1	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Streptococcus pyogenes	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 2	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Bordetella pertussis	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 3	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Mycobacterium tuberculosis	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 4	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Legionella pneumophila	2.0 x 10 <sup>6</sup> IFU/mL
Enterovirus EV71	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Mycoplasma pneumoniae	2.0 x 10 <sup>6</sup> U/mL
Respiratory syncytial virus	9x10 <sup>5</sup> TCID <sub>50</sub> /mL	Haemophilus influenzae	2.0 x 10 <sup>6</sup> IFU/mL
Rhinovirus	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Candida albicans	2.0 x 10 <sup>6</sup> IFU/mL
Influenza A virus (H1N1)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Staphylococcus aureus	2.0 x 10 <sup>6</sup> IFU/mL
Influenza A virus (H3N2)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Pseudomonas aeruginosa	8.7x10 <sup>7</sup> CFU/mL
Influenza B virus (Yamagata)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Escherichia coli	2.0 x 10 <sup>6</sup> IFU/mL
Influenza B virus (Victoria)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	SARS-coronavirus	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Adenovirus C1 Ad. 71	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Pooled human nasal wash (Normal Flora)	

## Interference

2. Rapid SARS-CoV-2 Antigen Test Card has tested samples with common microorganism. The results showed that these microorganisms had no effect on the specificity of the assay up to the listed concentration.

Microorganisms	Concentrations	Microorganisms	Concentrations
Human coronavirus 229E	2.0x10 <sup>5</sup> CID <sub>50</sub> /mL	MERS-coronavirus	2.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL
Human coronavirus OC43	2.0x10 <sup>5</sup> TCID <sub>50</sub> /mL	Chlamydia pneumoniae	2.0 x 10 <sup>6</sup> IFU/mL
Human coronavirus NL63	2.0x10 <sup>5</sup> TCID <sub>50</sub> /mL	Streptococcus pneumoniae	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 1	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Streptococcus pyogenes	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 2	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Bordetella pertussis	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 3	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Mycobacterium tuberculosis	2.0 x 10 <sup>6</sup> IFU/mL
Parainfluenza virus 4	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Legionella pneumophila	2.0 x 10 <sup>6</sup> IFU/mL
Enterovirus EV71	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Mycoplasma pneumoniae	2.0 x 10 <sup>6</sup> U/mL
Respiratory syncytial virus	9x10 <sup>5</sup> TCID <sub>50</sub> /mL	Haemophilus influenzae	2.0 x 10 <sup>6</sup> IFU/mL
Rhinovirus	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Candida albicans	2.0 x 10 <sup>6</sup> IFU/mL
Influenza A virus (H1N1)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Staphylococcus aureus	2.0 x 10 <sup>6</sup> IFU/mL
Influenza A virus (H3N2)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Pseudomonas aeruginosa	8.7x10 <sup>7</sup> CFU/mL
Influenza B virus (Yamagata)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Escherichia coli	2.0 x 10 <sup>6</sup> IFU/mL
Influenza B virus (Victoria)	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	SARS-coronavirus	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Adenovirus C1 Ad. 71	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL	Pooled human nasal wash (Normal Flora)	

## ACCURACY

2019-nCoV Antigen Rapid Test Kit

	RT-PCR		
	Confirmed	Excluded	Total
Positive	53	4	57
Negative	2	200	202
Total	55	204	259

Sensitivity = 53/55 x 100% = 96.36% (95% CI: 91.32% - 99.99%)  
Specificity = 200/204 x 100% = 98.04% (95% CI: 96.14% - 99.94%)  
Accuracy = 253/259 x 100% = 97.68% (95% CI: 95.83% - 99.54%)

## Endogenous Substances

Rapid SARS-CoV-2 Antigen Test Card has tested samples with common endogenous substances. The results showed that these substances had no affect on the specificity of the assay up to the listed concentration.

Substances	Concentrations	Substances	Concentrations
Whole Blood	4% v/v	Homeopathic (Alkalol)	10% v/v
Mucin	0.5% w/v	Nasitil Nasal Drop (Phenylephrine)	15% v/v
Tobramycin	0.0004% w/v	Illiadin Merck - Oxymetazoline	15% v/v
Mebucain (Menthol)	0.15% mg/mL	Allergo Comod Nasal Spray (Cromolyn)	15% v/v
VISION PUMP SPREY (Benzocain)	0.15% w/v	Avamys (Fluticasone Propionate)	5% v/v
Mupirocin	10 mg/mL	Zicam	5% w/v
Enfluvir (Osetamivir Phosphate)	0.5% w/v	DISINOL	Phenol 2 g/10ml
		(Phenol 2 g/10ml, Chlorbutanol 3.5 g/ml)	Chlorbutanol 3.5 g/ml

## LIMITATIONS

1. The test is limited to the qualitative detection of SARS-CoV-2 viral antigen in nasopharyngeal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined by this assay.
2. Proper specimen collection is critical, and failure to follow the procedure may give inaccurate results. Improper specimen collection, storage or repeated freezing and thawing of specimens can lead to inaccurate results.
3. A negative test result may occur if the level of antigen in a specimen is below the limit of detection of the test.
4. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
5. Negative test results do not rule out other potential non-SARS-CoV-2 viral infections. Negative results should be confirmed by molecular diagnosis if COVID-19 disease is suspected.
6. Positive test results do not rule out co-infections with other pathogens.
7. Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
8. The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.

## REFERENCES

1. Wu C, Liu Y, Yang Y, Zhang P, Zhong W, Wang Y, et al. (February 2020). "Analysis of therapeutic targets for SARS-CoV-2 and discovery of potential drugs by computational methods". Acta Pharmaceutica Sinica B. doi:10.1016.

## SYMBOL INDEX

IVD	In Vitro Diagnostic Use	Manufacturing Date	Expiry Date
Σ	Tests per Kit	Manufacturing	Keep Dry
LOT	Batch Number	Do not reuse	
2-30°C	Store between 2-30°C	Keep away from Sunlight	
REP	Authorized Representative	REF	Catalog#
i	See Instruction for Use		

## Manufacturer:

**ONCOSEM**  
Onkolojik Sistemler San. Ve Tic. A.Ş.  
Mustafa Kemal Mah. 2125. Sok. 6A/8  
Çankaya-Ankara | TÜRKİYE | +90 (312) 219 50 33  
info@oncosem.com | oncosem.com

## Distributor Europe:

**Casada International GmbH**  
Obermeiers Feld 3 | 33104 Paderborn  
Germany | +49 (0) 5254 - 930 98 0  
info@casada.com | www.casada.com

Creation date  
19.11.2020